



Variant calling using GATK

Khalid Mahmood 2022

https://www.melbournebioinformatics.org.au/tutorials/tutorials/variant_calling_gatk1/variant_calling_gatk1/



Workshop overview

1. Objectives

2. Introduction to genetic variants

Why we study variants?

Types of genetic variants

3. GATK

Overview

Best practices pipeline for variant calling

4. Resources and tools

Resources

Databases

Tools

Help

5. Workshop



1. Objectives

- We aim to cover:
 - Perform QC of sequencing data
 - Align raw reads to reference sequences
 - Perform alignment metric and generating a QC report
 - Prepare alignment data for variant calling
 - Identify simple variants using GATK HaplotypeCaller
 - Visualise simple variant data (VCF files)
 - Perform basic variant filtering



2. Introduction to genetic variants

- There are approximately 3 billion base pairs in the human genome
- Humans share 99.5% of DNA with other humans
- A **variant** is a difference between similar genomes.
- Usually a difference between DNA sequences we are studying and a reference genome.
- To describe a variant we give the location (genomic coordinates) and genetic change.

e.g. chr2 9834 A→G



2. Introduction to genetic variants

There is high degree of similarity but the human genome is large ~3 billion nucleotides.


This results in approximately 4-5 million variants between any individual and the reference genome.

These, seemingly small number of variations likely explains a significant proportion of phenotypic diversity among humans.



2. Types of genetic variants

Types of genetic variants:




Focus of this workshop: Calling short germline variants





3. GATK overview



 Genome Analysis Toolkit (GATK): software package to analyze highthroughput sequencing data




3. GATK overview

- Download available from
- https://github.com/broadinstitute/gatk/releases
- Tutorial version: GATK 4.2.0.0 (September 2021)
- Current version: GATK 4.2.2.0
- Explore GATK website gatk.broadinstitute.org
 - Tool index provides tools usage instructions
 - Technical documentation provides details on for example Algorithms
 - Forum provides access to Q&As and community discussions




3. GATK Best practices pipeline






3.1 Pre-processing



- A sequencing experiment results in a large volume of sequencing reads
- Reads are not mapped to a reference
- Reads can contain errors and technical artifacts
- e.g. a molecule sequenced multiple times will result in duplicate reads
- We need to filter and prepare the reads and the alignment data – ready for variant calling



3.1 Map to reference




BWA-MEM

- bwa mem -M -t 4 -R
 "@RG\tID:SRR622461.7\tSM:NA12878\tLB:ERR194147\tPL:ILLUMINA"
 <reference> sample_1.fastq sample_2.fastq > alignment.sam
- -M: inserts a tag to the alignment if non-primary alignment (required by GATK)
- -R: read group
- -t: threads or number of cpus
- <reference>: path to reference genome in fasta format and the BWA index files



3.1 Map to reference




BWA-MEM

- bwa mem -M -t 4 -R
 "@RG\tID:SRR622461.7\tSM:NA12878\tLB:ERR194147\tPL:ILLUMINA"
 <reference> sample 1.fastq sample 2.fastq > alignment.sam
- -R: read group contains information such as the sample name, library and flow cell.
- Refers to a set of reads generated from a single sequencing run in particular machine

@RG ID:SRR622461.7 SM:NA12878 LB:ERR194147 PL:ILLUMINA



3.1 Map to reference




- Output is a SAM/BAM file.
- SAM file specifications: https://samtools.github.io/hts-specs/SAMv1.pdf

```
Header
@HD
      VN:1.5 SO:coordinate
      ID:SRR622461.7 SM:NA12878
                                  LB:ERR194147
@RG
                                               PL:ILLUMINA
@PG
      ID:bwa PN:bwa VN:0.7.17-r1188 CL:bwa mem -M -t 4 -R
                                                                               flags/
Alignment
                   position
                                CIGAR
                                                       read
                                                                               metadata
ERR194147.45
              163
                   chr18 6576006 99 101M = 6578028 317
                                                       CATTTCT.... <B<<BBBBB...
                                                                             NM:i:0 MD:Z:101...
                                                                    PHRED
read name
             flag
                             mapping
                                          mate
                                          information
                             quality
                                                                    quality
```




3.1 Mark duplicates



- Mark Duplicates
- Identify reads that are non-independent measurement of sequence fragment
 - Same template of DNA sampled multiple times
 - PCR duplicates
- High sequence identify
- Align to same reference position



3.1 Mark duplicates




Mark Duplicates

- gatk MarkDuplicates -I sample.bam -O sample.dedup.bam -M sample.dedup.metrics.txt
- Recommended to be performed on reads per library or lane
- SAM flags are used to mark reads as duplicates
- Downstream GATK tools depend on these flags to assess support for variants and alleles




3.1 Base recalibration




- Base recalibration
- gatk tools BaseRecalibrator and ApplyRecalibration
- Performed per-sample to detect and correct for patterns of systematic errors in base quality scores.
- Evidenced by calculating metrics based on known variant locations
- Important for building reliable evidence for downstream analysis.


3. GATK Best practices pipeline







Software




HaplotypeCaller

CombineGVCFs/
GenomicsDBImport


GenotypeGVCFs






HaplotypeCaller

• gatk --java-options "-Xmx4g" HaplotypeCaller -R <reference.fa> -I input.bam -O output.g.vcf.gz -ERC GVCF








CombineGVCFs

 gatk CombineGVCFs R <reference.fa> --variant sample1.g.vcf.gz --variant sample2.g.vcf.gz -O cohort.g.vcf.gz

• Combine per samples gVCF files (produced by HaplotypeCaller) into a multi-sample gVCF file.








GenotypeGVCFs

gatk --java-options "-Xmx4g" GenotypeGVCFs -R <reference.fa>
 -V cohort.q.vcf.qz -O output.vcf.qz

• Combine per samples gVCF files (produced by HaplotypeCaller) into a multi-sample gVCF file.







- Output is a VCF file
- VCF file specifications https://samtools.github.io/hts-specs/VCFv4.2.pdf

Header

##FILTER=<ID=PASS,Description="All filters passed">
##contig=<ID=1,length=249250621>
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths for the ref and alt alleles in the order listed">


Variant record

#CHROM	POS	ID	REF	ALT	QUAL	FILTER	INFO	FORMAT	Sample1
1	567376		G	Α	146.3	PASS	AC=1;DP=55	GT:AD:DP	0/1:30,25:55





3.3 Variant Refinement



Refinement

- Variant callers are sensitive
- The aim here is to identify potential false positives and apply filters to remove those less likely to be real variants. Strategies include:
- 1. Variant quality score recalibration (using known sites)
- 2. Hard filtering on quality criteria
- 3. Annotation features

Variant record

#CHROM	POS	ID	REF	ALT	QUAL	FILTER	INFO	FORMAT	Sample gatk
1	567376		G	Α	146.3	PASS	AC=1;DP=55	GT:AD:DP	0/1:30,23.33



4. Resources and tools

- GATK resources bundle: collection of files for GATK based analysis working with human sequencing data.
- ftp://gsapubftp-anonymous@ftp.broadinstitute.org/bundle/hg38

```
1000G_omni2.5.hg38.vcf.gz
1000G_phase1.snps.high_confidence.hg38.vcf.gz
Axiom_Exome_Plus.genotypes.all_populations.poly.hg38.vcf.gz
dbsnp_146.hg38.vcf.gz
hapmap_3.3_grch38_pop_stratified_af.vcf.gz
hapmap_3.3.hg38.vcf.gz
Homo_sapiens_assembly38.dict
Homo_sapiens_assembly38.fasta
Homo_sapiens_assembly38.fasta.gz
Mills_and_1000G_gold_standard.indels.hg38.vcf.gz
```

4. Resources and tools

- BWA-MEM index
- bwa index Homo_sapiens_assembly38.fasta

Homo_sapiens_assembly38.fasta.amb Homo_sapiens_assembly38.fasta.amb Homo_sapiens_assembly38.fasta.ann Homo_sapiens_assembly38.fasta.bwt Homo_sapiens_assembly38.fasta.pac Homo_sapiens_assembly38.fasta.sa





4. Resources and tools

Tools name	function
FastQC	QC tools for raw sequencing reads
MultiQC	QC report aggregator (generates an HTML report)
GATK	Set of tools for variant calling
Picard	A command line tool to analysis and manipulate sequencing files
Samtools	Suite of tools for interacting with mapped sequencing reads (SAM/BAM/CRAM format)
BCFtools	Suite of tools for interacting with variant data (VCF/BCF formats)



4. Help

Tool documentation


• GATK forum

Online resources (e.g. Biostar)

GitHub for technical issues/discussions







Topic

Table of contents

Author Information

Overview

Learning Objectives

Description

Requirements and preparation

melbournebioinformatics/...

Mode of Delivery

Byobu-screen

Tutorial setting

The Genome Analysis Toolkit (GATK)

How this tutorial works

Tutorial contents table

Section 1: Map raw mapped reads to reference genome

- 1. Preparation and data import
- 2. Align genome

Section 2: Prepare analysis ready reads

- Sort SAM/BAM
- Mark duplicate reads
- 3. Base quality recalibration

Section 3: Variant calling


- Apply HaplotypeCaller
- 2. Apply CombineGVCFs
- 3. Apply GenotypeGVCFs

Section 4: Filter and prepare

Tutorial navigation



Command and output blocks '#' comments - do not run



Interactive sections
Notes, hints, exercises



Table of contents

Author Information

Overview

Learning Objectives

Description

Requirements and preparation

Mode of Delivery

Byobu-screen

Tutorial setting

The Genome Analysis Toolkit (GATK)

How this tutorial works

Tutorial contents table

Section 1: Map raw mapped reads to reference genome

- 1. Preparation and data import
- 2. Align genome

Section 2: Prepare analysis ready reads

- 1. Sort SAM/BAM
- 2. Mark duplicate reads
- 3. Base quality recalibration

Section 3: Variant calling

- 1. Apply HaplotypeCaller
- 2. Apply CombineGVCFs
- 3. Apply GenotypeGVCFs

Section 4: Filter and prepare analysis ready variants

- Variant Quality Score
 Recalibration
- Additional filtering
- 3. Final analysis ready VCF file

Section 5: Exporting variant data and visualisation

- 1. VariantsToTable
- 2. HTML report

Introductory material
Tutorial delivery and some instructions

Workshop content:

- 5 sections
- Each section covers a stage in the variant calling pipeline
- Each section has a text explain the process and links to relevant material
- Sections have multiple steps. Mostly have an input and an expected output file.
- This is a pipeline: input to a step is the output from a previous step



Workshop computers

We will be conducting the workshop on virtual machines

 Hosted on the University of Melbourne Research Cloud service and the ARDC Nectar Research Cloud infrastructure.

 Infrastructure for development and setup of the workshops machines by Simon Gladman







Workshop computers


Each participant should have a username and a password

- Each participant will be assigned one of the following two VM machines:
 - 115.146.84.252
 - 115.146.84.226
- Configuration:



Log on to the VMs

 Open a terminal window and on the command prompt type and enter:





Useful Linux commands

- Autofill on command line: Tab key
- Abort command: Ctrl-c
- List contents of a directory: ls -1
- What's the path to my current directory: pwd
- Change directory: cd <path/to/destination>
- Create a directory: mkdir <directory name>
- Copy a file: cp <source file> <destination path/name>
- Remove a directory: rmdir <directory name>
- Remove a file: rm <file name>
- Rename/move a file (this is not copying a file): mv <source file> <destination file>

- Open a text editor: nano
- Print file content (small files): cat <file
 name>
- Print file content (quick view): less <file name>
- Print file content (quick view/first 10 lines of a file): head <file name>
- Print file content (quick view/last 10 lines of a file): tail <file name>
- curl or wget: download a file from a URL (you will see this in other QIIME2 tutorials)
- Documentation for a command line tool: try man <tool name> OR <command name> --help



Workshop data

- Primary data: paired-end sequencing reads from the chr20
 - chr20:2677705-6631126
- Whole genome sequencing data
 - Female
 - Utah resident (European ancestry)
 - 1000 genomes project (NA12878)
- Other data from
 - A male and female
 - Utah resident (European ancestry)
 - 1000 genomes project (NA12891 and NA12892)

Byrska-Bishop, Marta et al. "High coverage whole genome sequencing of the expanded 1000 Genomes Project cohort including 602 trios". *bioRxiv*. (2021).



Byobu-screen

- A terminal multiplexer or a tool to to create multiple 'windows' in a single screen
- Improves stability of terminal sessions when connected to a remove computer
- List screen sessions: byobu-screen -ls
- Start new session: byobu-screen -S workshop
- Detach from screen to original window: Ctrl-a-d
- More details:
- https://www.melbournebioinformatics.org.au/tutorials/tutorials/variant_calling_gatk1/variant_calling_gatk1/#byobu-screen



Workshop

https://www.melbournebioinformatics.org.au/tutorials/tutorials/variant_calling_gatk1/variant_calling_gatk1/

